

ENDOSCOPIC SKULL BASE AND PITUITARY APPROACHES

**A STEP-BY-STEP GUIDE FOR SURGICAL INSTRUCTION
AND CADAVERIC DISSECTION**



Theodore H. SCHWARTZ

**Professor of Neurosurgery, Otorhinolaryngology and Neurology
and Neuroscience, Co-Director, Institute for Minimally Invasive Skull
Base and Pituitary Surgery, Weill Cornell Medical College**

New York-Presbyterian Hospital, New York, U.S.A.

Vijay K. ANAND

**Clinical Professor of Otorhinolaryngology
Co-Director, Institute for Minimally Invasive Skull Base
and Pituitary Surgery, Weill Cornell Medical College**

New York-Presbyterian Hospital, New York, U.S.A.

**ENDOSCOPIC SKULL BASE
AND PITUITARY APPROACHES**
**A STEP-BY-STEP GUIDE FOR SURGICAL INSTRUCTION
AND CADAVERIC DISSECTION**

Dr. Theodore H. SCHWARTZ

Professor of Neurosurgery, Otorhinolaryngology and Neurology and Neuroscience,
Co-Director, Institute for Minimally Invasive Skull Base and Pituitary Surgery

Correspondence address:

Department of Neurosurgery
Weill Cornell Medical College
New York Presbyterian Hospital
525 East 68th St., Box #99
New York, N.Y. 10021, U.S.A.
Phone: (212) 746-5620
Fax: (212)746-8947
schwarh@med.cornell.edu

<http://www.endoscopicskullbasesurgery.com>
<http://www.weillcornell.org/tschwartz>

Dr. Vijay K. ANAND

Clinical Professor of Otorhinolaryngology
Co-Director, Institute for Minimally Invasive Skull Base and Pituitary Surgery
Department of Otorhinolaryngology

Correspondence address:

772 Park Ave.
New York, NY 10021, U.S.A.
Phone: (212) 452-3005
Fax: (212) 452-3660
<http://www.sinusitis-solutions.com>
Vijayanandmd@aol.com

Acknowledgements:

We would like to thank Yvonne Bernschneider and Susan Jaffy for assistance with this publication. We would also like to thank Hugh Thomas for the medical illustration. Finally, we are indebted to our families for the enduring support and our patients for their trust.

Special thanks to Jonathan Roth and Ameet Singh for help with anatomic cadaver dissection.

Medical Illustrator:

Hugh Thomas
212-254-5463
hugh.thomas@verizon.net

Please note:

Medical knowledge is constantly changing. As new research and clinical experience broaden our knowledge, changes in treatment and therapy may be required. The authors and editors of the material herein have consulted sources believed to be reliable in their efforts to provide information that is complete and in accordance with the standards accepted at the time of publication. However, in view of the possibility of human error by the authors, editors, or publisher of the work herein, or changes in medical knowledge, neither the authors, editors, publisher, nor any other party who has been involved in the preparation of this work, can guarantee that the information contained herein is in every respect accurate or complete, and they cannot be held responsible for any errors or omissions or for the results obtained from use of such information. The information contained within this brochure is intended for use by doctors and other health care professionals. This material is not intended for use as a basis for treatment decisions, and is not a substitute for professional consultation and/or use of peer-reviewed medical literature.

Some of the product names, patents, and registered designs referred to in this booklet are in fact registered trademarks or proprietary names even though specific reference to this fact is not always made in the text. Therefore, the appearance of a name without designation as proprietary is not to be construed as a representation by the publisher that it is in the public domain.

**Endoscopic Skull Base and Pituitary Approaches
A Step-By-Step Guide for Surgical Instruction
and Cadaveric Dissection****Dr. Theodore H. SCHWARTZ**

Professor of Neurosurgery, Otorhinolaryngology and Neurology and Neuroscience, Co-Director, Institute for Minimally Invasive Skull Base and Pituitary Surgery

Dr. Vijay K. ANAND

Clinical Professor of Otorhinolaryngology
Co-Director, Institute for Minimally Invasive Skull Base and Pituitary Surgery
Department of Otorhinolaryngology

Correspondence addresses:**Dr. Theodore H. Schwartz, M.D., F.A.C.S.**

Professor of Neurosurgery, Otorhinolaryngology and Neurology and Neuroscience, Co-Director, Institute for Minimally Invasive Skull Base and Pituitary Surgery
Department of Neurosurgery
Weill Cornell Medical College
New York Presbyterian Hospital
525 East 68th St., Box #99
New York, N.Y. 10021, U.S.A.
Phone: (212) 746-5620
Fax: (212) 746-8947
schwarh@med.cornell.edu
www.endoscopicskullbasesurgery.com
http://www.weillcornell.org/tschwartz

Dr. Vijay K. Anand, M.D., F.A.C.S.

Clinical Professor of Otorhinolaryngology
Co-Director, Institute for Minimally Invasive Skull Base and Pituitary Surgery
Department of Otorhinolaryngology
772 Park Ave.
New York, NY 100211, U.S.A.
Phone: (212) 452-3005
Fax: (212) 452-3660
http://www.sinusitis-solutions.com
Vijayanandmd@aol.com

© 2014 Published by **Endo:Press®**, Tuttlingen
ISBN 978-3-89756-817-4, Printed in Germany
P.O. Box, D-78503 Tuttlingen, Germany
Phone: +49 74 61/1 45 90
Fax: +49 74 61/708-529
E-mail: Endopress@t-online.de

Editions in languages other than English and German are in preparation. For up-to-date information, please contact **Endo:Press®**, publisher Tuttlingen, Germany, at the address shown above.

Printed by:

Straub Druck+Medien AG,
D-78713 Schramberg, Germany

08.14-0.6

All rights reserved. No part of this publication may be translated, reprinted or reproduced, transmitted in any form or by any means, electronic or mechanical, now known or hereafter invented, including photocopying and recording, or utilized in any information storage or retrieval system without the prior written permission of the copyright holder.

Table of Contents

Preface	7
Introduction	7
Endoscopic Skull Base Corridors, Approaches and Targets	8
Transsphenoidal Corridor	9
Transnasal Corridor	16
Transethmoidal Corridor	19
Transmaxillary Corridor	22
Closure	25
Bibliography	26



Dr. Theodore H. Schwartz received his undergraduate and medical degrees from Harvard University where he graduated Magna Cum Laude. After completing his residency and chief residency in Neurosurgery at The Neurological Institute of New York at Columbia-Presbyterian Medical Center, Dr. Schwartz completed advanced fellowship training at Yale-New Haven Medical Center in the surgical treatment of brain tumors and epilepsy. Dr. Schwartz specializes in image-guided minimally invasive surgical techniques such as stereotaxis, endoscopy and intraoperative MRI, and has received numerous awards and fellowships including the prestigious van Wagenen Fellowship, awarded by the American Association of Neurological Surgeons and the von Humboldt Fellowship, awarded by the German Government. Dr. Schwartz was recently awarded the “Gentle Giant Award” by the Pituitary Network and is on the editorial board for the *Journal of Neurosurgery* and *World Neurosurgery*. Dr. Schwartz is a Professor of Neurosurgery, Otorhinolaryngology and Neurology and Neuroscience at Weill Cornell Medical Center, New York Presbyterian Hospital, surgical director of the Comprehensive Epilepsy Center, as well as co-director of the Institute for Minimally Invasive Skull Base and Pituitary Surgery and co-director of Surgical neuro-oncology. His recent textbook publications include *Practical Endoscopic Skull Base Surgery* (Plural Publishing Inc.) and *Endoscopic Pituitary Surgery* (Thieme).

Dr. Vijay K. Anand is a world renowned endoscopic sinus surgeon who was instrumental in developing image guidance in endoscopic sinus surgery and anterior skull base surgery. He graduated from Madras Medical College in India with honors and has won many awards in the field of medicine. He trained at Manhattan Eye Ear and Throat Hospital in Otolaryngology and currently is the Director of the Rhinology lab and postgraduate education in Rhinology at Weill Medical College of Cornell University and New York Presbyterian Hospital in New York. The Rhinology lab which has been funded by numerous private foundations is dedicated to research in basic and clinical sciences. He was the President of the American Rhinologic Society in 1995 and has been a pioneer in the development of endoscopic sinus surgery and its extended applications. He has published widely in the field of Rhinology including the recently published textbook on *Practical Endoscopic Skull Base Surgery* (Plural Publishing Inc.). Dr. Anand has been the Course Director and has conducted more than 35 courses in Advanced Endoscopic Sinus surgery at the Weill Medical College of Cornell University. He is the recipient of the Outstanding Teacher Award in Rhinology from the American Rhinological Society. He is a sought after speaker in the field of Rhinology and is a Clinical Professor of Otolaryngology at the Weill Medical College of Cornell University in New York as well as co-director of the Institute for Minimally Invasive Skull Base and Pituitary Surgery.

Preface

Endoscopic approaches to the pituitary and skull base are quickly becoming a standard of care in neurosurgery and otolaryngology. The basis for a comprehensive understanding of the applications and limitations of these approaches is best acquired in the laboratory, performing cadaveric dissections. After teaching and participating in several endoscopic skull base dissection courses, we felt there was a need for a dissection manual that could be helpful to guide the surgeon through the

various approaches in a step-by-step fashion. For this reason we have tried to make this manual very simple and illustrative. In addition, since the cadaveric anatomy never quite perfectly simulates real intraoperative conditions, we have linked each step with an intraoperative photo as a demonstration. The purpose of the intraoperative photographs are to assist the surgeon in making the cognitive transition from the cadaver laboratory to the operating room.

Introduction

The skull base lies at the anatomic boundary between the fields of neurosurgery and otolaryngology. Surgery in this region has always been a challenge for both disciplines. The success of endoscopic techniques in the management of inflammatory sinus disease has led to the next step of applying the endoscope to the resection of tumors of the skull base.

This laboratory manual is intended as a guide for cadaveric dissection which will serve as an introduction to the surgical exercises. We find it useful to think about the endoscopic skull base approaches as a combination of three factors – ① a target, ② a skull base approach and ③ a nasal corridor. The first aspect of the surgical plan is the target. We have defined 15 separate targets. They are – ① anterior fossa, ② olfactory groove, ③ sella, ④ suprasellar cistern, ⑤ lateral sphenoid sinus, ⑥ medial cavernous sinus, ⑦ lateral cavernous sinus, ⑧ orbital apex, ⑨ pterygopalatine fossa, ⑩ Meckel's cave, ⑪ infratemporal fossa, ⑫ petrous apex, ⑬ upper third of clivus, ⑭ lower

two-thirds of clivus, and ⑮ odontoid/ventral craniovertebral junction. Some targets have one possible approach, whereas other targets have multiple approaches. The second aspect of the approach involves an understanding of the possible corridors through which one passes on the way to the target. There are four corridors that define the endonasal endoscopic approaches: ① transnasal, ② transsphenoidal, ③ transethmoidal, ④ transmaxillary. These corridors correspond to the nasal sinuses and can be combined to reach a variety of targets. The link between the nasal corridor and the surgical target is the approach.


The purpose of the manual is to describe the various nasal corridors that are presently available to reach intracranial targets through the endonasal skull base approaches and to take the participant through each of these approaches. The division of the endonasal skull base approaches into corridors, approaches and targets provides a framework for instruction.

Endoscopic Skull Base Corridors, Approaches and Targets

Corridor	Approach	Target
Transnasal	Transcribriform Transclival Transodontoid	Anterior Fossa/Olfactory Groove Lower 2/3 of Clivus Odontoid/Craniovertebral Junction
Transsphenoidal	Transsellar Transtuberculum/ Transplanum Transclival Transcavernous	Sella Suprasellar Cistern Upper 1/3 of Clivus Medial Cavernous Sinus
Transthmoidal	Transfovea Ethmoidalis Transorbital* Transsphenoidal	Anterior Fossa Orbital Apex Cavernous Sinus
Transmaxillary	Transpterygoidal* Transpterygoidal* Transpterygoidal* Transpterygoidal* Transpterygoidal* Transpterygoidal*	Pterygopalatine Fossa Infratemporal Fossa Petrous Apex Lateral Sphenoid Sinus Lateral Cavernous Sinus Meckel's Cave

+ The transthmoidal transorbital approach involves opening the anterior and lateral walls of the sphenoid sinus.

* The transmaxillary transpterygoidal approach involves opening the ethmoid and sphenoid sinuses as well.



Transsphenoidal Corridor


The transsphenoidal corridor offers an approach to the sphenoid sinus which is the gateway to the sella, the planum sphenoidale, the suprasellar cistern, the intrasphenoidal clivus and medial cavernous sinus. Although a unilateral approach is feasible when

removing small pituitary tumors, the bilateral approach is critical for more extensive skull base approaches and eases visualization and instrument manipulation during cadaveric dissection.



Step 1 Advance the endoscope through each nares and identify the septum (**S**) medially, the inferior, middle (**MT**) and superior turbinates laterally and the choana (**CH**) inferiorly (**Fig. 1**).







Step 2 Identify the ostia of the sphenoid sinuses bilaterally (**OS**). These are found 1.5 cm above the choana (**CH**), just below the superior turbinates (**ST**). Below the ostium is the sphenothmoidal recess (**Fig. 2**). If a nasoseptal flap is being harvested, this must be done at the beginning of the operation (see Closure, p. 25). If not, proceed to **Step 3**.

Step 3 Remove the posterior nasal septum adjacent to the sphenoid rostrum.


Step 4 Complete a submucous resection of the nasal septum and harvest the bone of the vomer or the perpendicular plate of the ethmoid bone. These specimens can be used to reconstruct the skull base defect.

Step 5 Enlarge the ostia (**OS**) bilaterally with a mushroom forceps or a drill. Care must be taken not to damage the sphenopalatine artery when opening the ostium inferolaterally (**Fig. 3**).


Step 6 Identify the sphenoid rostrum and drill the bone until the sphenoid sinus is opened widely. Remove the intersphenoidal septae and identify the sella (**S**) and the clivus (**C**) inferiorly. The carotid protuberances (**CP**) overlie the vertical segment of the carotid artery on each side of the clivus (**Fig. 4**).




Step 7 Identify the anatomy on the lateral wall of the sphenoid sinus, such as the optic protuberance (**OP**), the carotid protuberance (**CP**) and the opticocarotid recess (**OCR**) (**Fig. 5**).





Step 8 Identify the tuberculum sellae (**TS**) and planum sphenoidale (**PS**) above the sella (**S**) (**Fig. 6**).





Step 9 Transsellar Approach: Drill the anterior wall of the sella (**S**) to expose the dura (**D**) overlying the pituitary gland. A large opening from carotid protuberance to carotid protuberance is helpful in removing large adenomas. Open the dura in a cruciate fashion to expose the pituitary gland (**Fig. 7**).








Step 10 Transtuberculum, Transplanum Approach:



Drill through the tuberculum sellae and planum sphenoidale to expose the dura of the tuberculum sellae (TSD) and planum sphenoidale (PSD). Open the dura above and below the superior intercavernous sinus (SICS) and then cauterize and cut the SICS to expose the optic chiasm (OC), optic nerves (ON) and arachnoid (A) overlying the gyrus rectus (GR). The anterior communicating artery (ACA) can be found above the optic chiasm and the pituitary stalk is below, a common site for craniopharyngiomas (CR). Advancing the endoscope through the floor of the third ventricle will expose the roof of the third ventricle, choroid plexus (CP), fornices (F) and foramina of Monro (FM) (Figs. 8a–d).








Step 11 Tuberculum sellae and planum meningiomas often involve the diaphragma sellae (**DS**), which has to be resected with care taken to preserve the pituitary stalk (**PS**). Meningiomas, which invade the optic canals must be removed in their entirety. Furthermore, opening of the optic canals is important in achieving this goal. Complete removal of planum meningiomas often exposes the lamina terminalis (**LT**) as well as the A1, anterior communicating artery (**ACA**) and A2 branches of the anterior cerebral artery. The view behind the stalk reveals the basilar artery (**BA**), posterior cerebral arteries (**PCA**), posterior communicating artery (**PCOM**), and third cranial nerves (**III**) arising from the mesencephalon. The optic tract (**OT**) and mamillary bodies (**MB**) also come into view (**Figs. 9a–c**).







Step 12 Transclival Approach (superior 1/3 of clivus): Drill the bone of the clivus to expose the clival dura (CD) below the pituitary dura (PD). The bone over the cavernous sinus (CS) and internal cerebral artery (ICA) can also be removed to expose the ICA in order to reach pathology like chordomas that extend behind the CS and ICA. Opening of the dura, which is done in the shape of a capital “I” to avoid damaging the sixth nerve, exposes the basilar artery (BA), third nerve (III), posterior cerebral arteries (PCA) and superior cerebellar arteries (SCA) in the interpeduncular cistern (Figs. 10a–c).







Step 13 To extend the transsphenoidal, transclival approach more laterally and superiorly to obtain a view behind the pituitary gland and carotid artery, one can remove the posterior clinoid (**PC**) after removing the bone over the pituitary dura (**PD**) and cavernous sinus (**CS**). The bone can be removed over the vertical portion of carotid artery (**ICA**) to mobilize the ICA laterally. In addition, the pituitary gland (**P**) is exposed and pushed medially to expose the dura over the posterior clinoid and sellar diaphragm, which can be removed for exposure of the third nerve, ambient cistern, medial temporal lobe and tentorium. This approach is useful for removal of the lateral and superior extend of petroclival meningiomas (**T**) (**Figs. 11a, b**).

Step 14 Transcavernous Approach: The bone between the optic nerve and carotid artery, or medial opticocarotid recess can be removed to expose the superomedial aspect of the cavernous sinus (**CS**). This opening can be extended inferolaterally to expose the carotid siphon (**ICA**) in the medial cavernous sinus. However, the cavernous sinus is best explored later in the dissection after the ethmoids have been opened through a corridor lateral to the middle turbinate and once the medial pterygoid plate has been removed (**Figs. 12a, b**).





Transnasal Corridor

The transnasal corridor lies medial to the middle turbinate and lateral to the septum. The superior border is the cribriform plate, the inferior border is the palate and posteriorly the nasal corridor extends to the choana, nasopharynx, inferior clivus and odontoid.



Step 1 Remove the perpendicular plate of the ethmoid bone to expose the cribriform plates bilaterally.

Step 2 Identify the vertical (**V**) attachment of the middle turbinate (**MT**). The cribriform plate (**CP**) lies medial to the vertical attachment of the middle turbinate and lateral to the septum (which was previously removed). This is a common site for meningoceles (**MC**) (**Fig. 13**).



Step 3 Transcribriform Approach: Remove the mucosa underlying the cribriform plate and drill the plate until paper thin and remove with a curette to expose the dura. Open the dura to expose the olfactory nerves (**ON**) bilaterally in the olfactory grooves (**Fig. 14**).



Step 4 Identify the choana and approach the nasopharynx.



Step 5 Laterally displace the inferior turbinate bilaterally with a Goldman bar.

Step 6 Completely remove the vomer. Identify the mucosa of the nasopharynx over the clivus and odontoid and identify the Eustachian tubes (**ET**) bilaterally (**Fig. 15a**).

Step 7 Drill the floor of the sphenoid sinus and identify the vidian nerves laterally. The vidian nerves represent the supero-lateral limits of the transnasal, transclival approach.

Step 8 Transclival Approach: Elevate a flap of mucosa and fascia at the back of the nasopharynx. This is done in an inverted U-shaped incision to be replaced at the end of the operation. The lateral limits of this flap are the Eustachian tubes (**ET**). The flap of basopharyngeal fascia is retracted downwards into the oropharynx to expose the clivus (**C**) and if necessary, the odontoid (**Figs. 15b, c**).






Step 9 Drill the clivus (**C**) to expose the clival dura (**CD**) (**Fig. 16a**).

Step 10 Open the dura in the shape of a capital “I” and cauterize back the edges (**Fig. 16b**).

Step 11 Identify the ventral aspect of the pons (**P**) and medulla (**M**) as well as the vertebral arteries and basilar artery, which may be visible, depending on how the pathology has displaced the normal anatomy.

Step 12 Transodontoid Approach: Extend the basopharyngeal fascia opening downwards to expose the odontoid (**O**) and ring of C1. Transect the atlantooccipital membrane, longus capitis and longus colli muscles (**LC**) (**Fig. 17**).

Step 13 The anterior arch of C1 can be removed to expose the dens which can be removed with a high speed drill after separating it from the apical and alar ligaments to expose the craniovertebral junction.





Transethmoidal Corridor

The transethmoidal corridor lies lateral to the middle turbinate and is the corridor to the fovea ethmoidalis, orbital apex and lateral sphenoidal sinus. The transethmoid approach is also useful in fully exposing the maxillary sinus and the transpterygoid approach



Step 1 Mobilize the middle turbinate medially.




Step 2 Working lateral to the septum (S) and middle turbinate (MT), identify the uncinus process (U) and starting with an uncinectomy and infundibulotomy, identify the ethmoidal bulla (EB) (Figs. 18a, b).



Step 3 Complete the uncinectomy superiorly and expose frontal recess.

Step 4 Penetrate the ethmoidal bulla and complete the ethmoidectomy to expose the fovea ethmoidalis (FE) (Fig. 19).

Step 5 Identify the anterior ethmoidal artery (AE) at the frontal recess (FR) and transect it after satisfactory clipping of the vessel. It is easily identified at the junction of the lamina papyracea (LP) and the frontal recess (Fig. 20).

Step 6 Complete the dissection posteriorly, inferiorly and medially to avoid injury to the orbital apex. The middle turbinates can be completely removed to facilitate the exposure.







Step 7 Transfovea Ethmoidalis Approach: Once a total bilateral ethmoidectomy is performed and the nasal septum and vertical attachments of the middle turbinates have been removed the cribriform plate dura (**CPD**) and fovea ethmoidalis (**FE**) and are exposed between the lamina papyracea (**LP**) of the medial orbits. Removal of the FE and opening of the dura expose the frontal lobes (**FL**) which can be seen in this example of removal of an olfactory groove meningioma (**OGM**) (**Figs. 21a, b**).

Step 8 Transcavernous Approach: Open the anterior wall of the sphenoid sinus and remove bone over the sella to expose the pituitary dura (**PD**). Additional removal of bone overlying the carotid artery will expose the cavernous



sinus (**CS**) more directly than the transsphenoidal corridor which is medial to the middle turbinate. The use of intraoperative Doppler is often useful to localize the carotid artery within the cavernous sinus (**Fig. 22**).

Step 9 Transorbital Approach: Remove the lamina papyracea to expose the medial wall of the orbit and the periorbital fat (**OF**). Care must be taken not to damage the medial rectus muscle. The bone removal can extend back into the sphenoid sinus to expose the orbital apex. The medial orbital apex generally presents to the lateral wall of the sphenoid sinus, although in 12–25% a posteriorly located ethmoid air cell or “Onodi cell” will contain the medial orbital apex (**Fig. 23**).




Transmaxillary Corridor

The transmaxillary corridor is used to reach the maxilla, pterygopalatine fossa, lateral sphenoid sinus and cavernous sinus, Meckel's cave, infratemporal fossa and petrous apex.



Step 1 At this point, the nasal septum (**S**) and middle turbinates (**MT**) have been removed. A total ethmoidectomy has been performed. Identify the opening of the maxillary sinus ostium (**MO**) (**Fig. 24**).




Step 2 Enlarge the posterior ridge of the maxillary sinus ostium and dissect the anterior vertical process of the pterygoid plate (**PP**) of the palatine bone to expose the sphenopalatine artery (**SPA**) and vidian nerve. Transect the sphenopalatine artery (**Fig. 25**).

Step 3 Transpterygoidal Approach: Remove the posterior plate of the palatine bone (**PP**) and drill the lateral wall of the sphenoid sinus as well as the posterior wall of the maxillary sinus (**MS**) to expose the pterygopalatine fossa (**Fig. 26**).

Step 4 The contents of the pterygopalatine fossa can be identified such as the origin of the sphenopalatine artery (**SPA**) emerging above the crista ethmoidalis (**CE**) as well as the descending palatine artery (**DPA**), posterior superior alveolar artery (**PSAA**) maxillary nerve (**Max N**) and infraorbital artery (**IOA**) and nerve (**ION**) (**Fig. 27**).



Step 5 Exposure of the posterior wall of the pterygopalatine fossa (**PPF**) reveals the vidian cana (**VN**) which can be followed posteriorly toward the ICA and the foramen rotundum (**FR**) which can be followed towards Meckel's cave and the middle fossa (**Fig. 28**).





Step 6 Further drilling of the lateral wall of the sphenoid sinus and the medial pterygoid bone exposes the lateral sphenoid sinus (**LSS**) and Meckel's cave (**MC**) superiorly and the petrous apex inferiorly (**Fig. 29**).


Step 7 Exposure of the infratemporal fossa (**ITF**) requires removal of the inferior turbinate, drilling of the crista ethmoidalis (**CE**) to fully expose the posterior wall of the maxillary sinus (**MAX**) and then drilling of the pterygo-maxillary fissure. The posterior wall of the maxillary sinus (**MAX**) is often thinned by the pathology and is easily fractured to reach pathology in the ITF (**Figs. 30a, b**).






Closure

A watertight closure of the skull base at the conclusion of endoscopic skull base surgery is as critical as the approach and the resection to prevent post-operative CSF leak and meningitis.

Step 1 Harvesting a Nasoseptal Flap: The nasoseptal flap (**NSF**) must be harvested at the beginning of the operation. The superior and inferior cuts are made a few millimeters below and above the junction of the septum and cribriform plate and the hard palate respectively. A third vertical cut is made as anterior as possible. Care is taken to preserve the vascular pedicle and the sphenopalatine artery (**Fig. 31**).



Step 2 The Gasket-Seal Closure: A piece of fascia lata (**FL**) is harvested, which is approximately 1 cm larger in diameter than the defect in the skull base. This graft is placed over the defect in the skull and then countersunk with a piece of vomer or Medpor® implant material (**P**) (Porex Surgical Inc., Newnan, GA, USA) which provides a rigid buttress for the closure. The edges of the **FL** stick out circumferentially like cauliflower providing a watertight “gasket” seal. The **NSF** is placed over the gasket seal so that the edges of the **NSF** extend beyond the **FL** and lie on the skull base. All mucosa must be removed from behind the **NSF** to prevent mucocoele formation. The **NSF** is held in place with a final layer of Tisseel (Baxter) or Duraseal (Covidien) (**DS**) (**Figs. 32a–c**).



Bibliography

- ANAND VK, SCHWARTZ TH, HILTZIK DH, KACKER A: Endoscopic transsphenoidal pituitary surgery with real-time intraoperative magnetic resonance imaging. *American Journal of Rhinology Am J Rhinol*. 2006 Jul-Aug;20(4):401-5.
- ANAND VK, SCHWARTZ TH: *Practical Endoscopic Skull Base Surgery*. Plural Publishing. San Diego, CA. 2007.
- ANAND VK, SCHWARTZ TH, McCOUL ED: Extended endoscopic skull base techniques. *Op. Tech. Otol-Head Neck Surg*. 22;4:253, 2011.
- ATTIA M, KANDASAMY J, JAKIMOVSKY D, BEDROSIAN J, ALIM M, ANAND VK, SCHWARTZ TH: The importance and timing of optic canal exploration and decompression during endonasal resection of tuberculum sellae and planum sphenoidale meningiomas. *Operative Neurosurgery* 71(1 Suppl Operative):58-67, 2012.
- ATTIA M, PATEL KS, KANDASAMY J, STEIG PE, SPINELLI HM, RIINA HA, ANAND VK, SCHWARTZ TH: Combined cranionasal surgery for sphenoidal meningiomas invading the paranasal sinuses, pterygopalatine and infratemporal fossa. *World Neurosurgery* 80;6:e367-73, 2013.
- BANU MA, SHIN BJ, KIM JH, WOODWORTH GF, BEDROSIAN JC, ANAND VK, SCHWARTZ TH: Low-dose intrathecal fluorescein and etiology-based graft choice in endoscopic endonasal closure CSF leaks. *Clinical Neurology and Neurosurgery* 116:28-34, 2014.
- BANU MA, RATHMAN A, SOUWEIDANE MM, ANAND VK, GREENFIELD JP, SCHWARTZ TH: Corridor-based endonasal endoscopic skull base surgery for pediatric skull base surgery with detailed radioanatomic measurements. *Operative Neurosurgery* (in press).
- BANU MA, GUERRERO A, RATHMAN A, SOUWEIDANE MM, ANAND VK, HEIER L, SCHWARTZ TH, GREENFIELD JP: Impact of skull base development on endonasal endoscopic surgical corridors. *Clinical article. Journal of Neurosurgery Pediatrics* 13;2:155-69, 2014.
- BEDROSIAN J, ANAND VK, SCHWARTZ TH: The endoscopic endonasal approaches for cranial base CSF leaks and meningo-encephaloceles. *World Neurosurgery* (in press).
- BEDROSIAN JC, GARCIA-NAVARRO V, McCOUL ED, ANAND VK, SCHWARTZ TH: Endoscopic balloon dilatation as an adjunct to extended endoscopic approaches to the skull base. *J. Neurosurg*. 116(6):1215-8, 2012.
- BEDROSIAN JC, RAITHATHA R, McCOUL ED, AKSELROD O, ANAND VK, SCHWARTZ TH: A prospective study of post-operative symptoms in sinonasal quality-of-life following endoscopic skull base surgery: dissociations based on specific symptoms. *Int Forum Allergy Rhinol* 3(8):664-9, 2013.
- BROWN SM, ANAND VK, TABAEE A, SCHWARTZ TH: Role of perioperative antibiotics in endoscopic skull base surgery. *Laryngoscope* 117: 1528-1532, 2007.
- BROWN SM, TABAEE A, SINGH A, SCHWARTZ TH, ANAND VK: Three-dimensional endoscopic sinus surgery: feasibility and technical aspects. *Otolaryngol Head Neck Surg* 138: 400-402, 2008.
- CAPPABIANCA P, ALFIERI A, de DIVITIIS E: Endoscopic endonasal transsphenoidal approach to the sella: towards functional endoscopic pituitary surgery (FEPS). *Minim Invasive Neurosurg* 41:66-73, 1998.
- CAPPABIANCA P, CAVALLO LM, COLAO A, de DIVITIIS E: Surgical complications associated with the endoscopic endonasal transsphenoidal approach for pituitary adenomas. *J Neurosurg* 97:293-298, 2002.
- CARRAU RL, JHO HD, KO Y: Transnasal-transsphenoidal endoscopic surgery of the pituitary gland. *Laryngoscope* 106:914-918, 1996.
- CAVALLO LM, CAPPABIANCA P, MESSINA A, ESPOSITO F, STELLA L, de DIVITIIS E, TSCHABITSCHER M: The extended endoscopic endonasal approach to the clivus and cranio vertebral junction: anatomical study. *Childs Nerv Syst*. 2007 Apr 6.
- CAVALLO LM, MESSINA A, CAPPABIANCA P, ESPOSITO F, de DIVITIIS E, GARDNER P, TSCHABITSCHER M: Endoscopic endonasal surgery of the midline skull base: anatomical study and clinical considerations. *Neurosurg Focus*. 2005 Jul 15;19(1):E2.
- CAVALLO LM, MESSINA A, GARDNER P, ESPOSITO F, KASSAM AB, CAPPABIANCA P, de DIVITIIS E, TSCHABITSCHER M: Extended endoscopic endonasal approach to the pterygopalatine fossa: anatomical study and clinical considerations. *Neurosurg Focus*. 2005 Jul 15;19(1):E5.
- COBB W, MAKOSCH G, ANAND VK, SCHWARTZ TH: Endoscopic transsphenoidal, transclival resection of an enterogenous cyst located ventral to the brainstem. *Neurosurgery* 67 (ONS Suppl 2): ONS518-519, 2010.
- de DIVITIIS E, CAPPABIANCA P, CAVALLO LM: Endoscopic transsphenoidal approach: adaptability of the procedure to different sellar lesions. *Neurosurgery* 51:699-705; discussion 705-697, 2002.
- de DIVITIIS E, CAVALLO LM, CAPPABIANCA P, ESPOSITO F: Extended endoscopic endonasal transsphenoidal approach for the removal of suprasellar tumors: Part 2. *Neurosurgery*. 2007 Jan;60(1):46-58; discussion 58-9.
- FRANK G, PASQUINI E, DOGLIETTO F, MAZZATENTA D, SCIARRETTA V, FARNETI G, CALBUCCI F: The endoscopic extended transsphenoidal approach for craniopharyngiomas. *Neurosurgery*. 2006 Jul;59 (1 Suppl 1):ONS75-83; discussion ONS75-83.
- FRANK G, SCIARRETTA V, MAZZATENTA D, FARNETI G, MODUGNO GC, PASQUINI E: Transsphenoidal endoscopic approach in the treatment of Rathke's cleft cyst. *Neurosurgery* 56:124-128; discussion 129, 2005.
- FRASER JF, ALLEN B, ANAND VK, SCHWARTZ TH: Three-dimensional neurostereotaxy: subjective and objective comparison to 2D. *Minim Invasive Neurosurg* 52: 25-31, 2009.


26. FRASER JF, ANAND VK, SCHWARTZ TH: Endoscopic biopsy sampling of tophaceous gout of the odontoid process. Case report and review of the literature. *J Neurosurg Spine* 7: 61–64, 2007.
27. FRASER JF, MASS AY, BROWN S, ANAND VK, SCHWARTZ TH: Transnasal endoscopic resection of a cavernous sinus hemangioma: technical note and review of the literature. *Skull Base* 18: 309–315, 2008.
28. FRASER JF, NYQUIST GG, MOORE N, ANAND VK, SCHWARTZ TH: Endoscopic endonasal transclival resection of chordomas: operative technique, clinical outcome, and review of the literature. *Journal of Neurosurgery* 112: 1061-1069, 2010.
29. FRASER JF, NYQUIST GG, MOORE N, ANAND VK, SCHWARTZ TH: Endoscopic endonasal minimal access approach to the clivus: case series and technical nuances. *Neurosurgery* 66 (ONS Suppl 1):ONS150-158, 2010.
30. GARCIA-NAVARRO V, ANAND VK, SCHWARTZ TH: The “Gasket-Seal” closure for extended endonasal endoscopic skull base surgery. Long-term efficacy in a large case series. *World Neurosurgery* 80;5:563-8, 2013.
31. GREENFIELD JP, LENG LZ, CHAUDHRY U, BROWN S, ANAND VK, SOUWEIDANE MM, SCHWARTZ TH: Combined simultaneous endoscopic transsphenoidal and endoscopic transventricular resection of a giant pituitary macroadenoma. *Minim Invasive Neurosurg* 51: 306–309, 2008.
32. GREENFIELD JG, KACKER A, BROWN TABAEE A, ANAND VK, SCHWARTZ TH: Endoscopic, endonasal, transethmoidal, transcribriform, transfovea ethmoidalis approach to the anterior skull base and anterior cranial fossa. *Neurosurgery* 62;5:883-892, 2010.
33. HOFSTETTER CP, SINGH A, ANAND VK, KACKER A, SCHWARTZ TH: The endoscopic, endonasal, transmaxillary transpterygoid approach to the pterygopalatine fossa, infratemporal fossa, petrous apex, and the Meckel cave. *J Neurosurg* 2009. (Epub ahead of print).
34. HOFSTETTER CP, SINGH A, ANAND VK, SCHWARTZ TH: Endoscopic transpterygoid approach to the pterygopalatine fossa, petrous apex, lateral sphenoid sinus and Meckel’s cave. *Journal of Neurosurgery* 113;5:967-74, 2010.
35. HOFSTETTER C, MANNAH RM, MUBITA L, ANAND VK, DEHDASHTI A, SCHWARTZ TH: Endoscopic endonasal approach for growth-hormone secreting pituitary adenomas. *Neurosurgical Focus* 29;4:E6, 2010.
36. HOFSTETTER C, SHIN BJ, MUBITA L, HUANG C, ANAND VK, BOOCKVAR JA, SCHWARTZ TH: Endoscopic endonasal surgery for functional pituitary adenomas. *Neurosurgical Focus* 30;4:E10, 2011.
37. HOFSTETTER CP, NANSZKO M, MUBITA L, TSOURIS J, ANAND VK, SCHWARTZ TH: Volumetric classification for giant pituitary macroadenomas predicts outcome and morbidity of endoscopic endonasal transsphenoidal surgery *Pituitary* 15(3):450-63, 2011.
38. HOFSTETTER, ANAND VK, SCHWARTZ TH: Endoscopic resection of pituitary adenomas *Op. Tech.Otol-Head Neck Surg.* 22;3:206-214, 2011.
39. JAKIMOVSKI D, ATTIA M, BONCI G, HOFSTETTER C, TSOURIS AJ, ANAND VK, SCHWARTZ TH: Incidence and significance of intraoperative CSF leak in endoscopic pituitary surgery using intrathecal fluorescein. *World Neurosurgery* (in press).
40. JHO HD, CARRAU RL, KO Y, DALY MA: Endoscopic pituitary surgery: an early experience. *Surg Neurol* 47:213–222; discussion 222–213, 1997.
41. JHO HD, CARRAU RL: Endoscopic endonasal transsphenoidal surgery: experience with 50 patients. *J Neurosurg* 87:44–51, 1997.
42. JHO HD, HA HG: Endoscopic endonasal skull base surgery: Part 1 – The midline anterior fossa skull base. *Minim Invasive Neurosurg* 47:1–8, 2004.
43. JHO HD: Endoscopic transsphenoidal surgery. *J Neurooncol* 54:187–195, 2001.
44. KASLIWAL MK, ANAND VK, LAVI E, SCHWARTZ TH: Endoscopic management of a rare case of nasal glioma in Meckel’s case in an adult: report of a case. *Min Invas Neurosurg* 53:191-193, 2010.
45. KASSAM A, SNYDERMAN CH, MINTZ A, GARDNER P, CARRAU RL: Expanded endonasal approach: the rostrocaudal axis. Part I. Crista galli to the sella turcica. *Neurosurg Focus.* 2005 Jul 15;19(1):E3. Review.
46. KASSAM A, SNYDERMAN CH, MINTZ A, GARDNER P, CARRAU RL: Expanded endonasal approach: the rostrocaudal axis. Part II. Posterior clinoids to the foramen magnum. *Neurosurg Focus.* 2005 Jul 15;19(1):E4.
47. KASSAM AB, GARDNER P, SNYDERMAN C, MINTZ A, CARRAU R: Expanded endonasal approach: fully endoscopic, completely transnasal approach to the middle third of the clivus, petrous bone, middle cranial fossa, and infratemporal fossa. *Neurosurg Focus.* 2005 Jul 15;19(1):E6.
48. KHAN OH, RAITHATHA R, ANAND VK, SCHWARTZ TH: Choosing the best approach in 3D endoscopic skull base surgery. *Otolaryngologia* 63;137-143, 2013.
49. KOMOTAR RJ, STARKE RM, RAPER DMS, ANAND VK, SCHWARTZ TH: The endoscope-assisted ventral approach compared with open microscope-assisted surgery for clival chordomas: a systematic review *World Neurosurgery* 76;3-4:318-27, 2011.
50. KOMOTAR RJ, STARKE RM, RAPER DMS, ANAND VK, SCHWARTZ TH: Endoscopic endonasal versus open repair of anterior skull base meningiomas: a systematic meta-analysis of outcomes. *World Neurosurgery* 77;5/6:713-724, 2012.
51. KOMOTAR RJ, STARKE RM, RAPER DMS, ANAND VK, SCHWARTZ TH: Endoscopic versus transsphenoidal microscopic and transcranial resection of craniopharyngiomas: a systematic meta-analysis of outcomes. *World Neurosurgery* 77;2:329-41, 2012.
52. KOMOTAR RJ, STARKE RM, RAPER DMS, ANAND VK, SCHWARTZ TH: Endoscopic endonasal compared with microscopic transsphenoidal and open transcranial resection of giant pituitary adenomas: a systematic review of outcomes. *Pituitary* 15(2):150-9, 2012.

53. KOMOTAR RJ, STARKE RM, RAPER DMS, ANAND VK, SCHWARTZ TH: Endoscopic endoscopic skull base surgery: a comprehensive comparison with open transcranial approaches. *Brit. J of Neurosurg.* 26(5): 637-48, 2012.
54. KOMOTAR RJ, STARKE RM, RAPER DM, ANAND VK, SCHWARTZ TH: Endoscopic endonasal versus open repair of anterior skull base CSF leak, meningocele and encephalocele. A systematic review of outcomes. *Journal of Neurological Surgery Part A. Central European Neurosurgery* 74(4):239-50, 2013.
55. KOMOTAR RJ, STARKE RM, RAPER DMS, ANAND VK, SCHWARTZ TH: Endoscopic endonasal versus open craniofacial resection of esthesioneuroblastoma *World Neurosurgery* 80;1-2:148-59, 2013.
56. KOMOTAR RJ, STARKE RM, RAPER DMS, ANAND VK, SCHWARTZ TH: Endonasal versus transoral odontoid resection: systematic meta-analysis of outcomes. *Innovative Neurosurgery* 1;1:37-47, 2013.
57. KULWIN C, SCHWARTZ TH, COHEN-GADOL AA: Endoscopic extended transsphenoidal resection of tuberculum sella meningiomas: nuances of neurosurgical technique. *Neurosurgery Focus* 35;6:E6, 2013.
58. LAUFER I, ANAND VK, SCHWARTZ TH: Endoscopic, endonasal extended transsphenoidal, transplanum, transtuberulum approach for resection of suprasellar lesions. *J Neurosurg.* 2007 Mar;106(3):400-6.
59. LAUFER I, GREENFIELD JP, ANAND VK, HARTL R, SCHWARTZ TH: Endonasal endoscopic resection of the odontoid process in a nonachondroplastic dwarf with juvenile rheumatoid arthritis: feasibility of the approach and utility of the intraoperative Iso-C three-dimensional navigation. Case report. *J Neurosurg Spine* 8: 376-380, 2008.
60. LEE DLY, McCOUL ED, ANAND VK, SCHWARTZ TH: Endoscopic endonasal access to the jugular foramen: defining the surgical approach. *J Neurological Surgery-Part B* 73;352-341, 2012.
61. LENG LZ, ANAND VK, HARTL R, SCHWARTZ TH: Endonasal endoscopic resection of an os odontoideum to decompress the cervicomedullary junction: a minimal access surgical technique. *Spine (Phila Pa 1976)* 34: E139-143, 2009.
62. LENG LZ, BROWN S, ANAND VK, SCHWARTZ TH: "Gasket-seal" watertight closure in minimal-access endoscopic cranial base surgery. *Neurosurgery* 62: ONSE342-343; discussion ONSE343, 2008.
63. LENG LZ, GREENFIELD JP, SOUWEIDANE MM, ANAND VK, SCHWARTZ TH: Endonasal endoscopic resection of craniopharyngiomas. Analysis of outcome measures including extent of resection, CSF leak, return to productivity and body mass index. *Neurosurgery* 70;1: 110-123, 2011.
64. LENG LZ, ANAND VK, SCHWARTZ TH: Endoscopic resection of craniopharyngomas *Op. Tech.Otol-Head Neck Surg.* Surgery 22;3:215-222, 2011.
65. LIOTTA D, KACKER A, SCHWARTZ TH, ANAND VK: Endoscopic management of juvenile nasopharyngeal angiofibromas. *Op. Tech.Otol-Head Neck Surg.* 22;4:281-284, 2011.
66. MASCARENHAS L, MOSHEL YA, SZENTERMAI O, SALEK AA, LENG LZ, HOFSTETTER CP, PLACANTONAKIS DG, TSOURIS AJ, ANAND VK, SCHWARTZ TH: The transplanum transtuberulum approaches for supr asellar and sellar-suprasellar lesions. Avoidance of CSF leak and lessons learned. *World Neurosurgery.* (in press).
67. McCOUL ED, ANAND VK, BEDROSIAN JR, SCHWARTZ TH: Endoscopic skull base surgery and its impact on sinonasal-related quality of life. *Int Forum All Rhin* 2(2):174-81, 2012.
68. McCOUL ED, SCHWARTZ TH, ANAND VK: Vascularized reconstruction of endoscopic skull base defects *Op. Tech. Otol-Head Neck Surg.* 22;3:232-236, 2011.
69. McCOUL ED, CHOW S, LEE D, ANAND VK, SCHWARTZ TH: Endoscopic endonasal resection of ventral skull base keratinaceous cysts. *Int Forum All Rhin* 2(3):258-63, 2012.
70. McCOUL ED, SCHWARTZ TH, ANAND VK: Endoscopic approach to the infratemporal fossa *Op. Tech.Otol-Head Neck Surg.* 22;4:285-290, 2011.
71. McCOUL ED, ANAND VK, SCHWARTZ TH: Endoscopic approaches to the cavernous sinus *Op. Tech.Otol-Head Neck Surg.* 22;4:263-269, 2011.
72. McCOUL ED, ANAND VK, SCHWARTZ TH: Significant improvements in health-related quality of life after endoscopic anterior skull base surgery. A prospective study. *J. Neurosurgery* 117:498-506, 2012.
73. McCOUL ED, ANAND VK, SCHWARTZ TH: Long-term effectiveness of a reconstruction protocol using the nasoseptal flap after endoscopic skull base surgery. *World Neurosurgery* 81;1:136-143, 2014.
74. NYQUIST GG, ANAND VK, MEHARA S, KICKER A, SCHWARTZ TH: Endoscopic repair of anterior skull base non-traumatic cerebrospinal fluid leaks, meningoceles and encephaloceles. *J. Neurosurgery* 112:1070-1077, 2010.
75. NYQUIST GG, ANAND VK, SINGH A, TABAEE A, SCHWARTZ TH: The Janus flap. The bilateral nasoseptal flap for anterior skull base reconstruction. *Oto Head Neck Surg* 142:327-331, 2010.
76. NYQUIST GG, ANAND VK, BROWN S, SINGH A, TABAEE A, SCHWARTZ TH: Middle turbinate preservation in endoscopic transsphenoidal surgery of the anterior skull base. *Skull Base* 20;5:343-347, 2010.
77. NYQUIST GG, ANAND VK, SCHWARTZ TH: Endoscopic management of cerebrospinal fluid rhinorrhea *Op. Tech. Otol-Head Neck Surg.* 22;3:229-231, 2011.
78. PATEL KS, KOMOTAR R, SZENTIRMAI O, MOUSSAZADEH N, RAPER DM, STARKE RM, ANAND VK, SCHWARTZ TH: Case-specific protocol to reduce CSF leakage after endonasal endoscopic surgery. *J. Neurosurgery* 119;3:661-8.2013.

79. PLACANTONAKIS DG, TABAEE A, ANAND VK, HILTZIK D, SCHWARTZ TH: Safety of low-dose intrathecal fluorescein in endoscopic cranial base surgery. *Neurosurgery* 61: 161–165; discussion 165–166, 2007.
80. RAITHATHA R, McCOUL ED, WOODWORTH GF, SCHWARTZ TH, ANAND VK: Endoscopic endonasal approaches to the cavernous sinus. *Int Forum Allergy Rhinol.* 2(1):9-15, 2012.
81. RAJAPPA P, MARGETIS K, SIGOUNAS D, ANAND VK, SCHWARTZ TH, GREENFIELD JG: Case report: endoscopic endonasal trans-clival approach to a ventral pontine pediatric ependymoma. *J. Neurosurgery Peds* 12;5:465-8, 2013.
82. RAPER DMS, KOMOTAR RJ, STARKE RM, ANAND VK, SCHWARTZ TH: Endoscopic versus open approaches to the skull base: a comprehensive literature review. *Op. Tech.Otol-Head Neck Surg.* 22;4:302-307, 2011.
83. ROTH J, SINGH A, NYQUIST G, FRASER JF, BERNARDO A, ANAND VK, SCHWARTZ TH: Three-dimensional and 2-dimensional endoscopic exposure of midline cranial base targets using expanded endonasal and transcranial approaches. *Neurosurgery* 65: 1116–1128; discussion 1128–1130, 2009.
84. ROTH J, FRASER J, SINGH A, BENARDO A, ANAND VK, SCHWARTZ TH: Surgical approaches to the orbital apex: Comparison of endoscopic endonasal and transcranial approaches using a novel 3D endoscope. *Orbit* 30;1: 43-48, 2011.
85. SCHABERG MR, ANAND VK, SCHWARTZ TH, COBB W: Microscopic versus endoscopic transnasal pituitary surgery. *Curr Opin Otolaryngol Head Neck Surg.* 2010 Feb;18(1):8-14.
86. SCHABERG M, ANAND V, SCHWARTZ TH: 10 Pearls for safe endoscopic skull base surgery. *Oto Clin North Am* 43:945-953, 2010.
87. SCHWARTZ TH, ANAND VK: The endoscopic endonasal transsphenoidal approach to the suprasellar cistern. *Clin Neurosurg* 54: 226-235, 2007.
88. SCHWARTZ TH, FRASER JF, BROWN S, TABAEE A, KACKER A, ANAND VK: Endoscopic cranial base surgery: classification of operative approaches. *Neurosurgery* 62: 991–1002; discussion 1002–1005, 2008.
89. SCHWARTZ TH, STIEG PE, ANAND VK: Endoscopic transsphenoidal pituitary surgery with intra-operative magnetic resonance imaging. *Neurosurgery.*2006 Feb;58 (1 Suppl):ONS44–51; discussion ONS44–51.
90. SILVA D, ATTIA M, KANDASAMY J, ALIM M, ANAND VK, SCHWARTZ TH: Endoscopic endonasal posterior clinoidectomy. A technical note. *Surg. Neurol Internatl* 3:64, 2012.
91. SILVA D, ATTIA M, KANDASAMY J, ALIM M, ANAND VK, SCHWARTZ TH: The endoscopic endonasal transsphenoidal “above and below” approach to the retroinfundibular area and interpeduncular cistern – cadaveric study and case illustrations. *World Neurosurgery* (in press).
92. SINGH A, WESELL A, ANAND VK, SCHWARTZ TH: Surgical anatomy and physiology for the skull base surgeon *Op. Tech.Otol-Head Neck Surg.* 22;3:184-193, 2011.
93. SOLARI D, MAGRO F, CAPPABIANCA P, CAVALLO LM, SAMII A, ESPOSITO F, PATERNO V, de DIVITIIS E, SAMII M: Anatomical study of the pterygopalatine fossa using an endoscopic endonasal approach: spatial relations and distances between surgical landmarks. *J Neurosurg.* 2007 Jan;106(1):157–63.
94. TABAEE A, ANAND VK, BARRON Y, HILTZIK DH, BROWN SM, KACKER A, MAZUMDAR M, SCHWARTZ TH: Endoscopic pituitary surgery: a systematic review and meta-analysis. *J Neurosurg* 111: 545–554, 2009.
95. TABAEE A, ANAND VK, BARRON Y, HILTZIK DH, BROWN SM, KACKER A, MAZUMDAR M, SCHWARTZ TH: Predictors of short-term outcomes following endoscopic pituitary surgery. *Clin Neurol Neurosurg* 111: 119–122, 2009.
96. TABAEE A, ANAND VK, BROWN SM, LIN JW, SCHWARTZ TH: Algorithm for reconstruction after endoscopic pituitary and skull base surgery. *Laryngoscope.* 2007 Jul;117(7):1133–7.
97. TABAEE A, ANAND VK, CAPPABIANCA P, STAMM A, ESPOSITO F, SCHWARTZ TH: Endoscopic management of spontaneous encephalocele of the lateral sphenoid sinus. *Journal of Neurosurgery* 112;5:1070-7, 2010.
98. TABAEE A, ANAND VK, FRASER JF, BROWN SM, SINGH A, SCHWARTZ TH: Three-dimensional endoscopic pituitary surgery. *Neurosurgery* 64: 288–293; discussion 294–285, 2009.
99. TABAEE A, PLACANTONAKIS DG, SCHWARTZ TH, ANAND VK: Intrathecal fluorescein in endoscopic skull base surgery. *Otolaryngol Head Neck Surg.* 2007 Aug;137(2):316–20.
100. TABAEE A, NYQUIST G, ANAND VK, SINGH A, KACKER A, SCHWARTZ TH: Palliative endoscopic surgery in advanced sinonasal and anterior skull base neoplasms. *Otolaryngol Head Neck Surg.* 2010 Jan;142(1):126-8.
101. WOODWORTH GF, MCCOUL ED, GREENFIELD JG, ANAND VK, SCHWARTZ TH: Endoscopic management of anterior cranial fossa meningiomas *Op. Tech.Otol-Head Neck Surg.* 22;4:254-263, 2011.
102. WOODWORTH GF, PATEL KS, SHIN B, BURKHARDT J-K, TSOURIS AJ, MCCOUL ED, ANAND VK, SCHWARTZ TH: Surgical outcomes using a medial-to-lateral endonasal endoscopic approach to pituitary adenomas invading the cavernous sinus. *Journal of Neurosurgery* (in press).

HOPKINS® Telescopes

for use in Endoscopically-Assisted Micro Neurosurgery (EAM)
Diameter 4 mm, length 18 cm



28132 AA – FA



0°

28132 AA

HOPKINS® Straight Forward Telescope 0°,
enlarged view, diameter 4 mm, length 18 cm,
autoclavable,
fiber optic light transmission incorporated,
color code: green



30°

28132 BA

HOPKINS® Forward-Oblique Telescope 30°,
enlarged view, diameter 4 mm, length 18 cm,
autoclavable,
fiber optic light transmission incorporated,
color code: red



45°

28132 FA

HOPKINS® Forward-Oblique Telescope 45°,
enlarged view, diameter 4 mm, length 18 cm,
autoclavable,
fiber optic light transmission incorporated,
color code: black



28132 FVA/CVA



45°

28132 FVA

HOPKINS® Forward-Oblique-Telescope 45°,
enlarged view, diameter 4 mm, length 18 cm,
autoclavable,
connection for fiber optic light cable upwards,
fiber optic light transmission incorporated,
color code: black



70°

28132 CVA


HOPKINS® Telescope 70°,
enlarged view, diameter 4 mm, length 18 cm,
autoclavable,
connection for fiber optic light cable upward,
fiber optic light transmission incorporated,
color code: yellow

It is recommended to check the suitability of the product for the intended procedure prior to use.


HOPKINS® Telescopes

for use in Endoscopically-Assisted Micro Neurosurgery (EAM)

Diameter 2.7 mm, length 18 cm



Diameter 4 mm, length 30 cm



Endoscope Holder

for Fixation of Flexible and Rigid Endoscopes



28272 RKB


28272 RKB **Holding System**, autoclavable,
with quick release coupling KSLOCK,
including:

Rotation Socket
to clamp to the OR table,
for European and US standard rails,
with lateral clamp for height and
angle adjustment of the articulated stand

Articulated Stand,
reinforced version, L-shaped,
with one central clamp for all five joint functions,
height 48 cm, swivel range 52 cm,
with quick release coupling KSLOCK (female)

Clamping Jaw,
metal, with axial intake,
clamping range 4.8 up to 12.5 mm,
with quick release coupling KSLOCK (male),
for use with instrument and telescope sheaths

KARL STORZ CLEARVISION® II System for intra-operative irrigation of the telescope lens



)*Optional Accessories:

MTP 031229-10 Single-use tubing set.
For use with Karl Storz CLEARVISION® II, sterile, 10 per pack

















Submit your order to:


mtp medical technical promotion gmbh,
Take-Off Gewerbepark 46, D-78579 Neuhausen ob Eck

KARL STORZ CLEARVISION® II

Irrigation Sheath for use with CLEARVISION® II System

Irrigation Sheath, proximally reinforced for use with Adjustable Holder 28272 RKB				Compatible HOPKINS® Telescopes			
							
Detail	Order No.	Outer Diameter	Working length	Order No.	View	Outer Diameter	Working length
	7230 AS	4.8 x 6.0 mm	14 cm	7230 AA	0°	4.0 mm	18 cm
	7230 BS	4.8 x 6.0 mm	14 cm	7230 BA	30°	4.0 mm	18 cm
	7230 FS	4.8 x 6.0 mm	14 cm	7230 FA	45°	4.0 mm	18 cm
	7230 CS	4.8 x 6.0 mm	14 cm	7230 CA	70°	4.0 mm	18 cm
	7220 AS	3.7 x 4.8 mm	10 cm	7220 AA	0°	3.0 mm	14 cm
	7220 BS	3.7 x 4.8 mm	10 cm	7220 BA	30°	3.0 mm	14 cm
	7220 FS	3.7 x 4.8 mm	10 cm	7220 FA	45°	3.0 mm	14 cm
	7220 CS	3.7 x 4.8 mm	10 cm	7220 CA	70°	3.0 mm	14 cm
	7219 AS	3.5 x 4.7 mm	14 cm	7229 AA	0°	2.7 mm	18 cm
	7219 BS	3.5 x 4.7 mm	14 cm	7229 BA	30°	2.7 mm	18 cm
	7219 FS	3.5 x 4.7 mm	14 cm	7229 FA	45°	2.7 mm	18 cm
	7219 CS	3.5 x 4.7 mm	14 cm	7229 CA	70°	2.7 mm	18 cm

Probes, Elevators, Knives and Currettes



- 629820 **Probe**, double-ended, maxillary sinus ostium seeker, ball-shaped ends diameter 1.2 and 2 mm, length 19 cm
- 628712 KUHN-BOLGER **Frontal Sinus Curette**, 55° curved, oval, forward cutting, length 19 cm
- 479100 COTTLE **Elevator**, double-ended, semisharp and blunt, graduated, length 20 cm
- 474000 FREER **Elevator**, double-ended, semisharp and blunt, length 20 cm
- 474001 FREER **Suction Elevator**, with stylet, length 19 cm
- 629830 KUHN **Frontal Ostium Seeker**, double-ended, No. 6, both sides curved 77°, one tip straight, other tip reverse angle, length 22 cm
- 628714 KUHN-BOLGER **Frontal Sinus Curette**, 90° curved, oval, forward cutting, length 19 cm
- 628001 **Sickle Knife**, pointed, length 19 cm
- 28164 KK de DIVITIIS-CAPPABIANCA **Scalpel**, with retractable blade, length 23 cm including:
Handle
Outer Sheath
Micro Knife, sickle-shaped
- 28164 EC CASTELNUOVO **Elevator**, double-ended, blunt end angled, semisharp end slightly curved, graduated, length 26 cm




28164 KK

28164 EC

RHINOFORCE® II Ethmoid Forceps

working length 17 cm



649100 B



649100 B **BLAKESLEY RHINOFORCE® II Ethmoid Forceps**, straight, size 0, with cleaning connector, working length 16 cm



649101 B **Same**, size 1



649102 B **Same**, size 2



649110 B **BLAKESLEY-WILDE RHINOFORCE® II Ethmoid Forceps**, 45° curved upwards, size 0, with cleaning connector, working length 16 cm




649111 B **Same**, size 1



649123 B **TAKAHASHI RHINOFORCE® II Ethmoid Forceps**, spoon size 4 x 10 mm, with cleaning connector, working length 16 cm

GRÜNWALD-HENKE/CASTELNUOVO RHINOFORCE® II Nasal Cutting Forceps through-cutting



451000 B GRÜNWALD-HENKE **RHINOFORCE® II Nasal Cutting Forceps**, straight, through-cutting, tissue-sparing, BLAKESLEY shape, size 0, width 3 mm, with cleaning connector, working length 13 cm

451001 B **Same**, size 1, width 3.5 mm

451002 B **Same**, size 2, width 4 mm



451010 B CASTELNUOVO **RHINOFORCE® II Nasal Cutting Forceps**, end of sheath 25° upturned, through-cutting, with straight jaws, BLAKESLEY shape, width 3 mm, with cleaning connector, working length 13 cm





451500 B GRÜNWALD-HENKE **RHINOFORCE® II Nasal Cutting Forceps**, 45° upturned, through-cutting, tissue-sparing, BLAKESLEY shape, size 0, width 3 mm, with cleaning connector, working length 13 cm

451501 B **Same**, size 1, width 3.5 mm

451502 B **Same**, size 2, width 4 mm




451510 B CASTELNUOVO **RHINOFORCE® II Nasal Cutting Forceps**, end of sheath 25° upturned, through-cutting, jaws 45° upturned, BLAKESLEY shape, width 3 mm, with cleaning connector, working length 13 cm

STRUYCKEN RHINOFORCE® II Nasal Cutting Forceps**TAKAHASHI RHINOFORCE® Nasal Forceps**

SilCut® Instruments

Special features:

- Tactile instrument feedback
- Uniform patented force transmission
- Powerful resection under precise control
- Accurate incision due to small tolerances
- Special cutting geometry to prevent tissue from slipping
- Large aperture angle
- Flat jaws
- Through-cutting and backward-cutting versions also available



451021

GRÜNWALD-HENKE **SilCut® Nasal Cutting Forceps**, straight, through-cutting, extremely powerful resection, patented uniform force transmission for gently controlled cutting, new ergonomic handle design, BLAKESLEY shape, size 1, with cleaning connector, working length 13 cm



451521

Same, 45° upturned

RHINOFORCE® II Miniature Nasal Forceps

flat jaws, through-cutting



452831

RHINOFORCE® II Miniature Nasal Forceps,
with extra fine flat jaws, through-cutting, tissue-sparing,
straight sheath, straight jaws, width of cut 1.5 mm,
with cleaning connector, working length 13 cm



452832

Same, jaws 45° upturned



452833


Same, sheath curved 30°, straight jaws



452834

Same, sheath curved 30°, jaws 45° upturned

STAMMBERGER Antrum Punch sidebiting downward and forward cutting



459051


STAMMBERGER Antrum Punch,
right side downward and forward cutting,
working length 10 cm



459052

Same, left side downward
and forward cutting

STAMMBERGER Antrum Punch backward cutting



459010

STAMMBERGER RHINOFORCE® II Antrum Punch,
upside backward cutting, with cleaning connector,
working length 10 cm



459011


Same, right side backward cutting



459012

Same, left side backward cutting


STAMMBERGER Antrum Punch backward cutting



459016

STAMMBERGER **Antrum Punch**, backward cutting, sheath 360° rotatable, with fixing screw, dismantling, working length 10 cm, for use with Cleaning Adaptor 459015 LL

STAMMBERGER Antrum Punch pediatric size, backward cutting




459036


STAMMBERGER **Antrum Punch**, small pediatric size, slender, backward cutting, sheath 360° rotating, with fixing screw, dismantling, working length 10 cm, for use with Cleaning Adaptor 459015 LL

STAMMBERGER Circular Cutting Punch


Special features:


- **Unique design**
 - **For enlarging openings in the sphenoid frontal wall**
 - **Circular cutting punch mechanism allows cutting in a full circle of 360° without rotating the instrument as required with a conventional punch forceps**
 - **No interference with other instruments simultaneously used in the nose (e. g. endoscope, suction tube)**
 - **Available in 2 sizes: diameter 3.5 and 4.5 mm, punch head 4-fold LASER-welded**
 - **Integrated irrigation channel**
- **Multipurpose use:**
 - In addition to abrading procedures applied to the sphenoid frontal wall, bony ethmoid septa, pieces of nasal concha and other thin bony bridges can also be cut away
 - Extremely useful for treatment of choanal atresia
 - Blunt punching head reduces injuries
 - If used correctly by punching exclusively in the sagittal axis, traumata of vital structures, e. g. cranium, arteria carotis and optic nerve osseous canal are virtually impossible










 651055 STAMMBERGER **Punch**, circular cutting, for sphenoid, ethmoid and choanal atresia, diameter 3.5 mm, with cleaning connector, working length 18 cm

 651050 **Same**, diameter 4.5 mm

 651060 STAMMBERGER **Punch**, circular cutting, 65° upturned, for frontal sinus recess, diameter 3.5 mm, with cleaning connector, working length 17 cm


 651065 **Same**, diameter 4.5 mm

STAMMBERGER Punches**Egg-shaped tip, for opening the ethmoid cells and the sphenoid sinus**

	651057	STAMMBERGER Punch , egg-shaped tip, circular cut, 60° cutting direction from distal above to proximal below, tip diameter 3.5 mm, straight sheath, for sphenoid, ethmoid and choanal atresia, with cleaning connector, working length 18 cm
	651058	Same , circular cut 120°
	651053	STAMMBERGER Punch , egg-shaped tip, circular cut, 120° cutting direction from distal below to proximal above, tip diameter 4.5 mm, straight sheath, for sphenoid, ethmoid and choanal atresia, with cleaning connector, working length 18 cm
	651052	Same , circular cut, 60° cutting direction, tip diameter 4.5 mm
	651061	STAMMBERGER Punch , egg-shaped tip, circular cut, 90° cutting direction, tip diameter 3.5 mm, sheath 65° upturned, for frontal sinus recess, with cleaning connector, working length 17 cm
	651066	Same , tip diameter 4.5 mm

Frontal Sinus Punches

with link chain sheath, backward cutting



651521 -
651523



651521

Frontal Sinus Punch, with link chain sheath 70° upturned, backward cutting, to reduce the spina nasalis superior, small, jaws 2.5 x 2 mm, working length 13 cm




651522

Same, medium (standard size), jaws 3.5 x 3 mm



651523

Same, large, jaws 5.5 x 5 mm

KERRISON **Bone Punches**

662101 – 62112


- | | | | |
|---|--|--------|--|
| □ | | 662101 | KERRISON Bone Punch , detachable, rigid, 90° upbiting, not through-cutting, size 1 mm, working length 17 cm |
| □ | | 662102 | Same , size 2 mm |
| □ | | 662112 | KERRISON Bone Punch , detachable, rigid, 90° downbiting, not through-cutting, size 2 mm, working length 17 cm |



662120

- | | | | |
|---|--|--------|---|
| □ | | 662120 | KERRISON Bone Punch , detachable, rigid, upbiting 40° forward, size 0.7 mm, working length 17 cm |
|---|--|--------|---|

Punches



28164 MKA – 28164 MKB




28164 MKA **KERRISON Bone Punch**, detachable, rigid, upbiting 60° forward, size 1 mm, working length 17 cm



28164 MKB **Same**, size 2 mm

Suction Tubes




10383 B – BL

10380 C – E

- 10383 B **Suction Tube**, with cut-off hole, diameter 3 mm, working length 35 cm
- 10383 BL **Same**, diameter 5.5 mm
- 10380 C **Suction Tube**, diameter 2 mm, working length 25 cm
- 10383 D **Same**, diameter 3 mm
- 10383 E **Same**, diameter 4 mm


Suction Tubes and Antrum Cannulas



- 529207 **FRAZIER Suction Tube**, with cut-off hole and stylet, angled, outer diameter 7 Fr./2 mm, working length 10 cm, total length 17.5 cm
- 586031 v. **EICKEN Antrum Cannula**, LUER-Lock, with cut-off hole, long curved, outer diameter 3 mm, length 12.5 cm
- 586226 v. **EICKEN Antrum Cannula**, LUER-Lock, with cut-off hole, short curved, outer diameter 2.5 mm, length 12.5 cm
- 586241 **Same**, outer diameter 4 mm


- 641625 **Suction Tube**, for frontal sinus, with cut-off hole, LUER, outer diameter 2.5 mm, length 14.5 cm
- 662885 **FRANK-PASQUINI Suction Tube**, angular, tip curved upwards, ball end, with grip plate and cut-off hole, LUER, diameter 3 mm, working length 13 cm
- 662886 **Same**, tip curved downwards

Suction Tubes and Irrigation Tube



- 206600 **FISCH Suction and Irrigation Tube**, cylindrical, suction tube outer diameter 2.5 mm, irrigation tube outer diameter 2 mm, working length 9.5 cm
- 649183 **FERGUSON Suction Tube**, with cut-off hole and stylet, LUER, 10 Fr., working length 15 cm
- 663818 **Suction Tube**, angular, malleable, with round handle and cut-off hole, diameter 2 mm, working length 13 cm
- 649179 B **Suction Tube**, malleable, with elongated cut-off hole and stylet, LUER, 4 Fr., working length 15 cm

- 649180 B **Suction Tube**, malleable, with elongated cut-off hole and stylet, LUER, 6 Fr., working length 15 cm
- 649182 B **Suction Tube**, malleable, with elongated cut-off hole and stylet, LUER, 8 Fr., working length 15 cm
- 649183 B **Suction Tube**, malleable, with elongated cut-off hole and stylet, LUER, 10 Fr., working length 15 cm

TAKE-APART® Bipolar Forceps

28164 BDL **TAKE-APART® Bipolar Forceps**, with fine jaws, width 1 mm, distally angled 45°, vertical closing, outer diameter 3.4 mm, working length 20 cm, including:
Bipolar Ring Handle
Outer Sheath
Inner Sheath
Forceps Insert



28164 BDM **TAKE-APART® Bipolar Forceps**, with fine jaws, width 1 mm, distally angled 45°, horizontal closing, outer diameter 3.4 mm, working length 20 cm, including:
Bipolar Ring Handle
Outer Sheath
Inner Sheath
Forceps Insert



28164 BDD **TAKE-APART® Bipolar Forceps**, width 2 mm, distally angled 45°, horizontal closing, outer diameter 3.4 mm, working length 20 cm, including:
Bipolar Ring Handle
Outer Sheath
Inner Sheath
Forceps Insert



28164 BDG **TAKE-APART® TAN Bipolar Coagulation Forceps**, size 3.4 mm, working length 20 cm, including:
Bipolar Ring Handle
Outer Sheath
Inner Sheath
Forceps Insert



26176 LA **Bipolar High Frequency Cord**, with 2x 4 mm banana plug for **KARL STORZ** Coagulator 26020 XA/XB and Valleylab, length 300 cm

STAMMBERGER Bipolar Suction Forceps



461010



461010

STAMMBERGER **Bipolar Suction Forceps**, 15° upturned, with suction channel, for bipolar coagulation in paranasal areas, working length 12.5 cm, for use with Bipolar High Frequency Cord 847002 E or 847002 A/M/V/U



461015



461015








STAMMBERGER **Bipolar Suction Forceps**, 45° upturned, with suction channel, for bipolar coagulation in paranasal areas, working length 12.5 cm, for use with Bipolar High Frequency Cord 847002 E or 847002 A/M/V/U

High Frequency Cords

for use with STAMMBERGER Bipolar Suction Forceps Accessories



Bipolar High Frequency Cords

KARL STORZ Instruments	High Frequency Electrosurgery Units		
		847002 E	<p>Bipolar High Frequency Cord, for KARL STORZ Coagulator 26021 B/C/D, 860021 B/C/D, 27810 B/C/D, 28810 B/C/D, AUTOCON® system (50, 200, 350), AUTOCON® II 400 SCB system (111, 113, 115) and Erbe coagulator, T and ICC series, with two 2 mm cable sockets for KARL STORZ Bipolar Suction Forceps 461010, 461015 and Bipolar Forceps 8615 A/AS, 28164 BGK, length 450 cm</p>
		847002 M	<p>Bipolar High Frequency Cord, for Martin and Berchtold coagulator, with two 2 mm cable sockets for KARL STORZ Bipolar Suction Forceps 461010, 461015 and Bipolar Forceps 8615 A/AS, 28164 BGK, length 450 cm</p>
		847002 A	<p>Bipolar High Frequency Cord, with 2 x 4 mm banana plug for KARL STORZ coagulator 26020 XA/XB, with two 2 mm cable sockets for KARL STORZ Bipolar Suction Forceps 461010, 461015 and Bipolar Forceps 8615 A/AS, 28164 BGK, length 450 cm</p>
		847002 V	<p>Bipolar High Frequency Cord, for KARL STORZ AUTOCON® II 400 SCB system (112, 114, 116), Valleylab coagulator, with two 2 mm cable sockets for KARL STORZ Bipolar Suction Forceps 461010, 461015 and Bipolar Forceps 8615 A/AS, 28164 BGK, length 450 cm</p>
KARL STORZ Instruments	Standard Forceps Bipolar Cords		
		847002 U	<p>Bipolar Universal High Frequency Cord, one side with two 2 mm cable sockets for KARL STORZ Bipolar Suction Forceps 461010, 461015 and Bipolar Forceps 8615 A/AS, 28164 BGK, other side with standard pin for connection to all current forceps bipolar cords, length 40 cm</p>

MONTGOMERY-YOUNGS RHINOFORCE® II Clip Applicator for endonasal endoscopic sphenopalatine artery ligature

Winner of the **2006**
Cutler Surgical Prize



The **Cutlers' Surgical Prize** is one of the most prestigious annual awards for innovation in the design or application of surgical instruments or techniques. The prize as well as the Clarke medal are awarded by the surgical association and the president of the Royal College of Surgeons of England.



452650 A MONTGOMERY-YOUNGS RHINOFORCE® II Clip Applicator, for endonasal endoscopic sphenopalatine artery ligature, with suction channel, handle with spring, straight, with cleaning connector, working length 13 cm, for use with Titanium Clips 8665 T



452650 C **Same**, jaws angled to the right



452650 D **Same**, jaws angled to the left

Coagulation Ball Electrode





28164 ED



28164 ED **Coagulation Ball Electrode**, diameter 2 mm, laterally curved, working length 13 cm




28164 EF **Same**, diameter 4 mm

Fine and Delicate Instruments**Forceps, working length 15 cm****Forceps, working length 18 cm**

Fine and Delicate Instruments

Miniature Forceps, working length 15 cm



662251

662251

Miniature Forceps, straight, through-cutting, with fine flat jaws, width of cut 1 mm, working length 15 cm

662255

Same, curved to right

662256

Same, curved to left


662257

Same, curved upwards

662271

Grasping Forceps, straight, fine-serrated, working length 15 cm

Miniature Forceps, working length 18 cm



28164 GS

28164 GS

Miniature Forceps, straight, through-cutting, with fine flat jaws, bite 1 mm, working length 18 cm

28164 GR

Same, curved right

28164 GL


Same, curved to left

28164 GU

Same, curved upwards

663271

Grasping Forceps, straight, fine-serrated, working length 18 cm

Fine and Delicate Instruments**Scissors, working length 15 cm****Scissors, working length 18 cm**

Sellar Stage

Curettes, round spoon



28164 KA



28164 KA **Curette**, round spoon, tip slightly angled, size 1 mm, with round handle, length 25 cm



28164 KB **Curette**, round spoon, tip slightly angled, size 2 mm, with round handle, length 25 cm



28164 KC **Curette**, round spoon, tip slightly angled, size 3 mm, with round handle, length 25 cm



28164 KF **Curette**, round spoon, tip highly angled, size 2 mm, with round handle, length 25 cm



28164 KG **Curette**, round spoon, tip highly angled, size 3 mm, with round handle, length 25 cm

de DIVITIIS-CAPPABIANCA Suction Curettes, with styplet, round wire – basket-shaped



28164 RSB



28164 RSB CAPPABIANCA-de DIVITIIS **Suction Curette**, blunt, inner diameter 5 mm, tip angled 45°, LUER, length 25 cm

28164 RSC **Same**, inner diameter 7 mm



28164 RT CAPPABIANCA-de DIVITIIS **Suction Curette**, with basket, round, size 5 mm, rotatable tube, LUER, length 25 cm



28164 RU **Same**, size 6.5 mm

CASTELNUOVO Elevator



28164 EA



28164 EA CASTELNUOVO **Elevator**, double-ended, semisharp and blunt, length 26 cm











Sellar Stage

Curettes



28164 RP






Inner diameter
in mm:

	28164 RN	CAPPABIANCA-de DIVITIIS Ring Curette , with round wire, inner diameter 3 mm, tip angled 45°, with round handle, length 25 cm
	28164 RO	CAPPABIANCA-de DIVITIIS Ring Curette , with round wire, inner diameter 5 mm, tip angled 45°, with round handle, length 25 cm
	28164 RP	CAPPABIANCA-de DIVITIIS Ring Curette , with round wire, inner diameter 7 mm, tip angled 45°, with round handle, length 25 cm
	28164 RG	CAPPABIANCA-de DIVITIIS Ring Curette , with round wire, inner diameter 5 mm, tip angled 90°, with round handle, length 25 cm
	28164 RH	Same , inner diameter 7 mm
	28164 RB	CAPPABIANCA-de DIVITIIS Ring Curette , with round wire, inner diameter 3 mm, laterally curved sheath end, with round handle, length 25 cm
	28164 RA	Same , inner diameter 5 mm
	28164 RC	Same , inner diameter 7 mm
	28164 RD	CAPPABIANCA-de DIVITIIS Ring Curette , with round wire, inner diameter 5 mm, laterally curved 90° sheath end, with round handle, length 25 cm


Delicate Dissectors and Elevators



28164 DA

	28164 DA	Dissector , sharp, tip angled 45°, round spatula, with round handle, size 2 mm, length 25 cm
	28164 DB	Dissector , sharp, tip angled 45°, round spatula, with round handle, size 3 mm, length 25 cm
	28164 DF	Dissector , sharp, tip angled 15°, flat long spatula, with round handle, size 1.5 mm, length 25 cm
	28164 DS	Dissector , sharp, tip angled 15°, with round handle, size 2 mm, length 25 cm
	28164 DM	Dissector , sharp, straight tip, slightly curved spatula, with round handle, size 3 mm, length 25 cm

Antrum Curette, Frontal Sinus Curette, Antrum Cannula and Suction Tube



628701 **Antrum Curette**, round, length 19 cm
 628702 **Same**, oblong, small size
 628703 **Same**, large size
 628712 KUHN-BOLGER **Frontal Sinus Curette**,
 55° curved, oval, forward cutting,
 length 19 cm
 586025 v. EICKEN **Antrum Cannula**, LUER-Lock,
 long curved, outer diameter 2.5 mm,
 length 12.5 cm
 586030 **Same**, outer diameter 3 mm
 586040 **Same**, outer diameter 4 mm

586125 v. EICKEN **Antrum Cannula**, LUER-Lock,
 long curved, malleable, serrated grip plate,
 outer diameter 2.5 mm, length 12.5 cm
 586130 **Same**, outer diameter 3 mm
 529309 FRAZIER **Suction Tube**, with mandrel and
 cut-off hole, with distance marking at 5–9 cm,
 9 Fr., working length 10 cm

SEPEHRNIA Neurosurgical Micro-Instruments**Needle Holder and Forceps**

28164 NBC

28164 NBC **Micro Needle Holder**, bayonet-shaped,
jaws curved to left, 1 x 6 mm,
working length 10 cm

28164 PBE

28164 PBB **Micro Forceps**, bayonet-shaped, spoon,
2 mm, working length 10 cm28164 PBE **Same**, 4 mm spoon

28164 PBG

28164 PBG **Micro Forceps**, bayonet-shaped,
spoon horizontal, 2 mm, working length 10 cm28164 PBH **Same**, 4 mm spoon horizontal

SEPEHRNIA Neurosurgical Micro-Instruments

Scissors



28164 SBC



28164 SBA **Micro Scissors**, bayonet-shaped, sharp/sharp, cutting edges straight, working length 10 cm



28164 SBB **Same**, bayonet-shaped, sharp/sharp, left curved



28164 SBC **Same**, bayonet shaped, blunt/blunt, jaw straight



28164 SBD **Same**, bayonet shaped, sharp/sharp, jaw curved to right



28164 SBE **Same**, bayonet shaped, sharp/sharp, jaws horizontal

Dissectores and Ring Curettes

GAAB Recommended Instruments



28164 GFO



- 28164 GBO **Dissector**, bayonet-shaped, sharp, round spatula, tip angled upwards 45°, with round handle, size 3 mm, working length 15 cm
28164 GBU **Same**, tip angled downwards 45°



- 28164 GFO **Dissector**, sharp, flat long spatula, tip angled upwards 15°, with round handle, size 1.5 mm, working length 15 cm
28164 GFU **Same**, tip angled downwards 15°



28164 GGO



- 28164 GGO **Ring Curette**, bayonet-shaped, round wire, inner diameter 5 mm, tip angled upwards 90°, with round handle, working length 15 cm
28164 GGU **Same**, tip angled downwards 90°



- NEW** 28164 GKO **Ring Curette**, bayonet-shaped, blunt, tip angled upwards 45°, outer diameter 4 mm, working length 15 cm
NEW 28164 GKU **Same**, tip angled downwards 45°



- NEW** 28164 GLL **Ring Curette**, bayonet-shaped, blunt, tip angled to left 90°, outer diameter 3.3 mm, with round handle, working length 15 cm
NEW 28164 GLR **Same**, tip angled to right 90°

Micro Vascular Knife and Dissector



28164 GM

28164 GM **Micro Vascular Knife**, bayonet-shaped, curved downwards, length 18.5 cm



28164 DL

28164 DL **Dissector**, bayonet-shaped, sharp, curved to left, length 11 cm

28164 DR **Same**, curved to right

IMAGE1 S Camera System ^{NEW}



Economical and future-proof

- Modular concept for flexible, rigid and 3D endoscopy as well as new technologies
- Forward and backward compatibility with video endoscopes and FULL HD camera heads


- Sustainable investment
- Compatible with all light sources




Innovative Design

- Dashboard: Complete overview with intuitive menu guidance
- Live menu: User-friendly and customizable
- Intelligent icons: Graphic representation changes when settings of connected devices or the entire system are adjusted


- Automatic light source control
- Side-by-side view: Parallel display of standard image and the Visualization mode
- Multiple source control: IMAGE1 S allows the simultaneous display, processing and documentation of image information from two connected image sources, e.g., for hybrid operations




Dashboard



Live menu



Intelligent icons



Side-by-side view: Parallel display of standard image and Visualization mode


IMAGE1 S Camera System ^{NEW}




Brilliant Imaging

- Clear and razor-sharp endoscopic images in FULL HD
- Natural color rendition


- Reflection is minimized
- Multiple IMAGE1 S technologies for homogeneous illumination, contrast enhancement and color shifting




FULL HD image




CLARA



FULL HD image



CHROMA



FULL HD image



SPECTRA A*



FULL HD image



SPECTRA B**

* SPECTRA A: Not for sale in the U.S.

** SPECTRA B: Not for sale in the U.S.

IMAGE1 S Camera System ^{NEW}**IMAGE1 S**

TC 200EN

- TC 200EN* **IMAGE1 S CONNECT**, connect module, for use with up to 3 link modules, resolution 1920 x 1080 pixels, with integrated KARL STORZ-SCB and digital Image Processing Module, power supply 100–120 VAC/200–240 VAC, 50/60 Hz including:
- Mains Cord**, length 300 cm
 - DVI-D Connecting Cable**, length 300 cm
 - SCB Connecting Cable**, length 100 cm
 - USB Flash Drive**, 32 GB, USB silicone keyboard, with touchpad, US
- * Available in the following languages: DE, ES, FR, IT, PT, RU

Specifications:

HD video outputs	- 2x DVI-D - 1x 3G-SDI
Format signal outputs	1920 x 1080p, 50/60 Hz
LINK video inputs	3x
USB interface	4x USB, (2x front, 2x rear)
SCB interface	2x 6-pin mini-DIN

Power supply	100–120 VAC/200–240 VAC
Power frequency	50/60 Hz
Protection class	I, CF-Defib
Dimensions w x h x d	305 x 54 x 320 mm
Weight	2.1 kg

**For use with IMAGE1 S
IMAGE1 S CONNECT Module TC 200EN**

TC 300

- TC 300 **IMAGE1 S H3-LINK**, link module, for use with IMAGE1 FULL HD three-chip camera heads, power supply 100–120 VAC/200–240 VAC, 50/60 Hz, **for use with IMAGE1 S CONNECT TC 200EN** including:
- Mains Cord**, length 300 cm
 - Link Cable**, length 20 cm

Specifications:

Camera System	TC 300 (H3-Link)
Supported camera heads/video endoscopes	TH 100, TH 101, TH 102, TH 103, TH 104, TH 106 (fully compatible with IMAGE1 S) 22220055-3, 22220056-3, 22220053-3, 22220060-3, 22220061-3, 22220054-3, 22220085-3 (compatible without IMAGE1 S technologies CLARA, CHROMA, SPECTRA*)
LINK video outputs	1x
Power supply	100–120 VAC/200–240 VAC
Power frequency	50/60 Hz
Protection class	I, CF-Defib
Dimensions w x h x d	305 x 54 x 320 mm
Weight	1.86 kg

* SPECTRA A: Not for sale in the U.S.

** SPECTRA B: Not for sale in the U.S.

IMAGE1 S Camera Heads ^{NEW}



For use with IMAGE1 S Camera System
IMAGE1 S CONNECT Module TC 200EN, IMAGE1 S H3-LINK Module TC 300
and with all IMAGE1 HUB™ HD Camera Control Units



TH 100

TH 100 **IMAGE1 S H3-Z Three-Chip FULL HD Camera Head**, 50/60 Hz, IMAGE1 S compatible, progressive scan, soakable, gas- and plasma-sterilizable, with integrated Parfocal Zoom Lens, focal length $f = 15-31$ mm (2x), 2 freely programmable camera head buttons, for use with IMAGE1 S and IMAGE1 HUB™ HD/HD

Specifications:

IMAGE1 FULL HD Camera Heads	IMAGE1 S H3-Z
Product no.	TH 100
Image sensor	3x 1/8" CCD chip
Dimensions w x h x d	39 x 49 x 114 mm
Weight	270 g
Optical interface	integrated Parfocal Zoom Lens, $f = 15-31$ mm (2x)
Min. sensitivity	F 1.4/1.17 Lux
Grip mechanism	standard eyepiece adaptor
Cable	non-detachable
Cable length	300 cm



TH 104

TH 104 **IMAGE1 S H3-ZA Three-Chip FULL HD Camera Head**, 50/60 Hz, IMAGE1 S compatible, **autoclavable**, progressive scan, soakable, gas- and plasma-sterilizable, with integrated Parfocal Zoom Lens, focal length $f = 15-31$ mm (2x), 2 freely programmable camera head buttons, for use with IMAGE1 S and IMAGE1 HUB™ HD/HD

Specifications:

IMAGE1 FULL HD Camera Heads	IMAGE1 S H3-ZA
Product no.	TH 104
Image sensor	3x 1/8" CCD chip
Dimensions w x h x d	39 x 49 x 100 mm
Weight	299 g
Optical interface	integrated Parfocal Zoom Lens, $f = 15-31$ mm (2x)
Min. sensitivity	F 1.4/1.17 Lux
Grip mechanism	standard eyepiece adaptor
Cable	non-detachable
Cable length	300 cm

Monitors



9619 NB

9619 NB

19" HD Monitor,
color systems **PAL/NTSC**, max. screen
resolution 1280 x 1024, image format 4:3,
power supply 100–240 VAC, 50/60 Hz,
wall-mounted with VESA 100 adaption,
including:
External 24 VDC Power Supply
Mains Cord



9826 NB

9826 NB

26" FULL HD Monitor,
wall-mounted with VESA 100 adaption,
color systems **PAL/NTSC**,
max. screen resolution 1920 x 1080,
image format 16:9,
power supply 100–240 VAC, 50/60 Hz
including:
External 24 VDC Power Supply
Mains Cord

Monitors

KARL STORZ HD and FULL HD Monitors	19"	26"
Wall-mounted with VESA 100 adaption	9619 NB	9826 NB
Inputs:		
DVI-D	●	●
Fibre Optic	-	-
3G-SDI	-	●
RGBS (VGA)	●	●
S-Video	●	●
Composite/FBAS	●	●
Outputs:		
DVI-D	●	●
S-Video	●	-
Composite/FBAS	●	●
RGBS (VGA)	●	-
3G-SDI	-	●
Signal Format Display:		
4:3	●	●
5:4	●	●
16:9	●	●
Picture-in-Picture	●	●
PAL/NTSC compatible	●	●

Optional accessories:

9826 SF **Pedestal**, for monitor 9826 NB
9626 SF **Pedestal**, for monitor 9619 NB

Specifications:

KARL STORZ HD and FULL HD Monitors	19"	26"
Desktop with pedestal	optional	optional
Product no.	9619 NB	9826 NB
Brightness	200 cd/m ² (typ)	500 cd/m ² (typ)
Max. viewing angle	178° vertical	178° vertical
Pixel distance	0.29 mm	0.3 mm
Reaction time	5 ms	8 ms
Contrast ratio	700:1	1400:1
Mount	100 mm VESA	100 mm VESA
Weight	7.6 kg	7.7 kg
Rated power	28 W	72 W
Operating conditions	0–40°C	5–35°C
Storage	-20–60°C	-20–60°C
Rel. humidity	max. 85%	max. 85%
Dimensions w x h x d	469.5 x 416 x 75.5 mm	643 x 396 x 87 mm
Power supply	100–240 VAC	100–240 VAC
Certified to	EN 60601-1, protection class IPX0	EN 60601-1, UL 60601-1, MDD93/42/EEC, protection class IPX2

Accessories for Video Documentation



- 495 NL **Fiber Optic Light Cable,**
with straight connector, diameter 3.5 mm,
length 180 cm
- 495 NA **Same,** length 230 cm

Cold Light Fountain XENON 300 SCB



- 20133101-1 **Cold Light Fountain XENON 300 SCB**
with built-in antifog air-pump, and integrated
KARL STORZ Communication Bus System SCB
power supply:
100–125 VAC/220–240 VAC, 50/60 Hz
including:
Mains Cord
Silicone Tubing Set, autoclavable, length 250 cm
SCB Connecting Cable, length 100 cm
- 20133027 **Spare Lamp Module XENON**
with heat sink, 300 watt, 15 volt
- 20133028 **XENON Spare Lamp,** only,
300 watt, 15 volt

Cold Light Fountain XENON NOVA® 300



- 20134001 **Cold Light Fountain XENON NOVA® 300,**
power supply:
100–125 VCA/220–240 VAC, 50/60 Hz
including:
Mains Cord
- 20132028 **XENON Spare Lamp,** only,
300 watt, 15 volt

KARL STORZ AIDA® compact NEO advanced

Brilliance in documentation



AIDA compact NEO:
Recording screen



AIDA compact NEO:
Patient data



AIDA compact NEO:
Review screen

Data Acquisition

Still images, video sequences and audio comments can easily be recorded during an examination or intervention by pressing the on-screen button, activating the footswitch, or pressing the camera head button.

All captured data are displayed on the right-hand side as a thumbnail preview to ensure the data have been generated. Patient data can be entered via an onscreen or standard keyboard. The system also offers the possibility to transfer all relevant patient data via a DICOM worklist or a link to the hospital information system (HIS) without requiring manual entry in the patient entry screen.

Flexible Review, Data Storage and Efficient Data Export

Captured still images or video files can easily be viewed, edited, or deleted on-screen before final storage. KARL STORZ AIDA® compact NEO efficiently stores all recorded data on DVD, CD, USB stick, external/internal drive, the relevant network and/or on a FTP server. It is also possible to save the data directly on the PACS and/or HIS servers via HL7/DICOM. Data that cannot be stored successfully remains in a cache until final archiving is possible.

Special Features:

- **SD and HD signal support:**
 - Y/C (S-Video)
 - Composite input
 - DVI-D input
- **Picture-in-Picture function:**
Display of channel 2 (SD) in channel 1 (FULL HD)
- **Resolution:**
 - Still images 1920 x 1080 and SD
 - Videos 1080p, 720p and SD
- Interface package (DICOM/H7) included
- NEO Secure security software
- Recommended applications:
 - Universal (cart or OR1™ installation)



20 0409 13-EN* **KARL STORZ AIDA® compact NEO advanced**

Documentation system for digital storage of still images, video sequences and audio files, power supply 115/230 VAC, 50/60 Hz

* Available in the following languages:
DE, ES, FR, IT, PT, PL, RU, DK, SE, JP, CN

Equipment Cart



UG 220

UG 220

Equipment Cart

wide, high, rides on 4 antistatic dual wheels equipped with locking brakes 3 shelves, mains switch on top cover, central beam with integrated electrical subdistributors with 12 sockets, holder for power supplies, potential earth connectors and cable winding on the outside,

Dimensions:

Equipment cart: 830 x 1474 x 730 mm (w x h x d),

shelf: 630 x 510 mm (w x d),

caster diameter: 150 mm

including:

Base module equipment cart, wide

Cover equipment, equipment cart wide


Beam package equipment, equipment cart high

3x **Shelf,** wide

Drawer unit with lock, wide

2x **Equipment rail,** long

Camera holder



UG 540

UG 540

Monitor Swivel Arm,

height and side adjustable, can be turned to the left or the right side, swivel range 180°, overhang 780 mm, overhang from centre 1170 mm, load capacity max. 15 kg, with monitor fixation VESA 5/100, for usage with equipment carts UG xxx

Recommended Accessories for Equipment Cart



UG 310

UG 310

Isolation Transformer,
200 V–240 V; 2000 VA with 3 special mains socket,
expulsion fuses, 3 grounding plugs,
dimensions: 330 x 90 x 495 mm (w x h x d),
for usage with equipment carts UG xxx



UG 410

UG 410

Earth Leakage Monitor,
200 V–240 V, for mounting at equipment cart,
control panel dimensions: 44 x 80 x 29 mm (w x h x d),
for usage with isolation transformer UG 310



UG 510

UG 510

Monitor Holding Arm,
height adjustable, inclinable,
mountable on left or right,
turning radius approx. 320°, overhang 530 mm,
load capacity max. 15 kg,
monitor fixation VESA 75/100,
for usage with equipment carts UG xxx

Notes:

**WITH COMPLIMENTS OF
KARL STORZ—ENDOSKOPE**